

Advancing Biomarker Discovery: Characterization of the P2 Technology for Plasma Samples for Pharma and Biobank, Use and Application to NSCLC



Arthur Viodé¹, Sandra Schär¹, Polina Shichkova¹, Tobias Treiber¹, Christopher Below¹, Nicole Barlow¹, Jan Muntel¹, Wouter van Bergen¹, Angelica Rigutto², Alessandra Curioni-Fontecedro², Carina Sihlbom Wallem³, Annika Thorsell³, Stina Lassesson³, Johannes Fuchs³, Amaury Lachaud¹, Yuehan Feng¹, Lukas Reiter¹, Roland Bruderer¹

¹Biognosys AG, Schlieren, Switzerland; ²Faculty of Science and Medicine, University of Fribourg, Fribourg, Switzerland; ³Proteomics Core Facility at University of Gothenburg, part of BioMS and SciLifeLab, Gothenburg, Sweden

ABSTRACT PRESENTATION NUMBER: MP 141

INTRODUCTION

Plasma is the most widely collected biofluid and an invaluable source of biomarkers. The analysis of plasma using discovery mass spectrometry-based proteomics faces challenges as the 22 most abundant proteins constitute more than 99% of the total protein content. This is hindering the detection of lower abundant proteins, potentially disease-relevant biomarkers. To overcome this, we developed and optimized an enrichment workflow, termed P2 Plasma Enrichment System, based on protein corona formation. Protein corona formation leads to a reduction in dynamic range enabling the detection of lower abundant proteins. We characterized P2 in terms of pre-analytical variation, robustness and quantification. Finally, we applied it to find predictive biomarkers in a multicentric phase II clinical trial (SAKK17/18) in patients with non-small cell lung cancer (NSCLC).




Figure 1 P2 Enrichment System is based on protein corona formation. Using an optimized particle surface and optimal buffer conditions reproducible protein corona is formed. This workflow is single-well and high-throughput. Adapted from Gopinath, 2019, Sci.Rep.¹

METHODS




Figure 2 P2 enrichment workflow created with BioRender.com

We performed P2 plasma protein enrichment starting with 100µL of plasma followed by tryptic digestion. Samples were processed in a semi-automated fashion using a KingFisher Flex. We first characterized our workflow robustness over 21 plates (**Figure 3**), labs (**Figure 4**), batches of P2 materials (**Figure 5**) and finally to blood collection methods (**Figure 7**). Further, we prepared and analyzed 95 samples from phase II clinical trial (SAKK17/18) in patients with NSCLC (**Figure 8**). Plasma samples were collected at baseline (n=33) and at two post-dose time points (C1D8, n=31 and C2D1, n=31) (**Figure 9–11**). The samples were acquired in data-independent acquisition (DIA) either on a timsTOF HT (Bruker) or on an Exploris 480 mass spectrometer (Thermo Scientific) equipped with FAIMS Pro. Data analysis was performed using Spectronaut® 19 (Biognosys) using directDIA™.

RESULTS: WORKFLOW CHARACTERIZATION

PLASMA QC SHOWS GREAT STABILITY ACROSS 21 96-WELL PLATES AND 2.5 MONTHS




Figure 3 One plasma QC was added to each 96-well plate and analyzed using SN 19 directDIA (A) protein groups (PG) identified in each QC (plate) (B) Coefficient of variation (CV) % of protein groups across the 21 plates. **There is no drop in protein identified over 2.5 months and the identifications are stable.**

INTERLABORATORY REPRODUCIBILITY: BIOGNOSYS AG, BIOGNOSYS INC AND GOTHENBURG UNIVERSITY




Figure 4 Interlaboratory comparison using the same plasma QC. Samples were prepared independently in each laboratory. Lab 1 is Biognosys AG, Lab 2 is Biognosys Inc. and Lab 3 is the University of Gothenburg. (A) Protein groups identified by each laboratory (n=3) (B) CV% of protein groups. Note that for Lab 2 the QCs were processed and analyzed on different days. **The P2 workflow is transferable to other labs.**

BATCH TO BATCH REPRODUCIBILITY




Figure 5 Four batches of the P2 particles were tested in quadruplicate. All batches performed similarly in terms of protein groups identified and CV. **There is very low variability between P2 particles batches.**

LINEARITY AND QUANTITATION WAS EVALUATED USING A PROTEIN SPIKE IN EXPERIMENT




Figure 6 Protein were spiked in plasma 6 different amount 0.4ng, 2ng, 10ng, 50ng, 250ng and 1250ng, prepared in triplicates. (A) IL6 and (B) is the phosphorylated peptides of Tau p217. **In each experiment the correlation R² was >0.95, highlighting the linearity of the enrichment.**

P2 PERFORMANCE ACCROSS BLOOD COLLECTION TYPES




Figure 7 Five different plasma collection methods were evaluated: ACD (Acid Citrate Dextrose), CPD (Citrate Phosphate Dextrose), NaCitrate (Sodium Citrate), K2EDTA and K3EDTA. For each plasma type, samples from 4 match individuals were collected. **Blood collection is a key pre-analytical factor to control for.**

RESULTS: BIOMARKER DISCOVERY




Figure 8 SAKK 17/18 Origin phase II multicentric trial overview.

PLASMA PROTEIN CHANGES DIFFERENTIATE TREATMENT RESPONSE BEFORE AND AFTER DOSING




Figure 9 Protein differential abundance analysis (ROTS²) reveals significantly up- and downregulated plasma proteins comparing partial response (PR) and progressive disease (PD) at (A) pre-dose/baseline; (B) C1D8; (C) C2D1. (D) Principal component analysis demonstrating patient stratification based on clinical outcomes at baseline using the 16-protein signature panel.

WGCNA IDENTIFIED PLASMA PROTEINS CORRELATED WITH PROGRESSION-FREE SURVIVAL (PFS)




Figure 10 Using weighted gene co-expression network analysis (WGCNA³), PFS driver proteins were identified based on module membership and gene significance. (A) WGCNA cluster dendrogram; (B) PFS-correlated (cyan module) driver proteins.

COX SURVIVAL ANALYSIS




Figure 11 Cox survival analysis. For each protein we fit the multivariate model as Surv(os_time, os) ~ age + sex + protein abundance.

Using ML (**Figure 9**) and WGCNA (**Figure 10**) models we identified 23 predictive biomarkers. Using these biomarkers we performed Cox survival analysis Seven proteins (**Figure 11**) showed significant associations with patient survival: two upregulated in responders (green) and five in non-responders (red), some of which have been previously reported. Notably, CCL17, upregulated in responders, has been linked to enhanced immune cell infiltration and antitumor activity in the LUAD tumor microenvironment⁴.

REFERENCES

- Gopinath et al., Sci.Rep. 2019
- Suomi et al., Plos Comp Bio. 2017
- Langfelder and Horvath, BMC Bioinformatics. 2008
- Ye et al., Front Cell Dev Biol. 2022

CONCLUSIONS

- P2 workflow enables unbiased and deep coverage of proteins in the plasma proteome: Robustness was demonstrated over 21 plate and 3 independent laboratories Linearity is preserved P2 works best in EDTA plasma
- An integrated pipeline, starting from sample collection, sample analysis and bioinformatics has been established for a phase II NSCLC trial.
- A panel of circulating protein signatures has been identified to be closely associated with progression-free survival.
- The presented approach provides a powerful platform for identifying predictive biomarkers, enabling more precise patient stratification and potentially improving clinical outcomes.

CONTACT

For project consultation and contract research services, please contact:

Arthur Viodé, PhD
Senior Scientist I, R&D
E: arthur.viode@biognosys.com
E: services@biognosys.com

