



## CCI User Rules

All users are required to read and understand the following CCI User Rules and confirm their acceptance in order to provide convenient and reliable access to high-quality imaging for as many people as possible while protecting both the users and the instruments from the consequences of any misuse.

### 1. Acknowledgement of the CCI in Publications

Acknowledgements are of the utmost importance for the CCI, as publications from our users are a major hallmark for the impact and importance of a facility. They prove the effective use of the CCI resources to university executives and research networks, and are a prerequisite for the success of the CCI in evaluation processes and grant applications.

When you (your colleagues/collaborators) use data acquired, analysed or processed at the CCI in a publication, or use them in a grant application, the **CCI must be acknowledged as part of VR-supported National Microscopy Infrastructure**. For example: *"We acknowledge the Centre for Cellular Imaging at the University of Gothenburg and the National Microscopy Infrastructure, NMI (VR-RFI 2019-00217) for providing assistance in microscopy"*. Publications encompass, but are not limited to, scientific papers, Bachelor-, Master- or PhD-theses as well as contributions to conferences.

### 2. First time user

- Both you, as a User, and your Principal Investigator (PI) or Project Leader need to set up an account in the **CrossLab**, our core-management system. The **CrossLab** is used for all CCI resource booking, registration, usage tracking, service requests, reporting and billing.
- The **CCI Registration Form**, which is found in CrossLab under "Request Services", has to be filled and signed by both the user and his/her project leader before becoming a CCI user.
- All research projects of national rank, i.e. where the advanced equipment at CCI facility is used to answer complex research questions, must be registered in the web-portal of the **National Microscopy Infrastructure**.
- The user must complete a training by the CCI team before getting access to the microscopes and CCI preparation lab, without exceptions! Usage of the CCI preparation lab is subject to the rules outlined in **Correct Labwork Procedure at the CCI and Risk Assessment** documents (contact CCI staff).

**IMPORTANT: Users cannot train new-users, and untrained guests may not operate the instruments/use the preparation lab!**

### 3. Training

Based on their training and experience, users are divided into the following categories, and have different access rights to the CCI:

- **First level:** During training and assisted usage.
- **Second level:** Inexperienced trained users, who might need support and technical assistance. Access restricted to office hours (08.00-17.00).
- **Third level:** Experienced trained users. No access restrictions.

### 4. General rules

- **You are responsible** for the system during your booked session. **Inform the CCI Staff immediately if there are any problems with the equipment!** If you are also uncertain about something or cannot remember how to operate the system, please ask the CCI staff for advice.
- Outside working hours (08.00-17.00) and on weekends/Holidays, **only third level users** (see Training) are allowed in the facility (microscope use and sample preparation).
- No more than two people can work at the same workstation/microscope simultaneously.
- Avoid the microscope illumination sources (lasers, lamps, filaments, EHT) being turned on and off more than once a day. Always look at the booking calendar to see if someone will come after you, and **follow microscope-specific guidelines from the CCI staff regarding how to turn on/off illumination sources and the whole system.**
- **Keep the equipment and all working spaces clean!** After you finish, i.e. clean the objectives (LM) and remove all your belongings.
- If an instrument is not available due to technical issues or maintenance, it is blocked for booking. Affected users will be informed by e-mail or telephone and will be supported by CCI staff to perform their experiments on an alternative instrument if possible.
- Access to CCI preparation lab is subject to users being approved to profit the equipment during working hours and out of hours by the CCI staff.

### 5. Booking and pricing

- **Booking** of the equipment in CrossLab before use **is mandatory!**
- **Real-time usage** of all equipment connected to CCI computers is **tracked via CrossLab** and is matched against the resource schedule in CrossLab (*Your group will be charged for the maximum time!*). Upon arrival to the equipment, you will use a web-based tool called the "Kiosk" to track your equipment usage.
- **No-show** (makes a reservation and neither uses nor cancels it) without any explanation prior to the booking will be charged for the whole booking period.
- The microscopes will be charged per half hour, with a minimum of one hour, i.e. working from 9:30-10:00 will count as an hour.
- The systems can be reserved up to **2 weeks** in advance and the maximum length of a reservation is **4 hours** (in one session) per day, unless you have requested additional time from the CCI staff. **Booking 4 hours-slot for more than three days in the same week** need to be approved by the CCI staff. Long-term time-lapse experiments (more than 5 hours) should be scheduled after 18:00.
- Cancellations:
  - should be done as early as possible, minimum **24 hours** in advance **to avoid being charged.**
  - **on the same day, either due to sickness or experiment failure**, will be charged for 1hr.
  - Remember to **Log your usage**, please note that if the machine is booked and usage is not logged, it will be considered a no-show and charged for the whole booking period.
- To ensure that all users are up to date with our rules and with how to use and handle the microscopes, **your booking permission will be put on hold if you have not been active for more than a year.** If that happens, you will need to schedule a refresher session with the CCI staff.

## 6. Data handling

- While acquiring the images, save them in a folder titled with your full name in the **Users** folder (on the desktop of each microscope computer).
- At the end of the imaging session, copy your images immediately to the **CCI transfer server**. It is **strictly forbidden to use USB sticks!** Data transfer is allowed in only one direction: FROM the CCI computers TO user backup storage solution. The CCI transfer server is used only for transfer and short-term storage (max. 3 weeks), and must be regularly emptied after the images have been backed up.
- It is not allowed to use the CCI transfer server for long-term storage or upload anything, apart from images or data acquired at the CCI.
- Once the data have been transferred from the CCI server to your storage solution **you have full responsibility for the data**. This includes both storage and backup as well as correct handling (i.e. that no image manipulation/misconduct is performed) and interpretation.

**IMPORTANT: All user images will be deleted on the microscope and analysis computers at the end of the month and only users that acquired data within the last 5 working days will be contacted!**

## 7. Safety Policies

- **Human tissue:** If working with human material, then you need to present a copy of the ethical permit before bringing any samples to the facility.
- **Animal Use:** Users planning to use live animals in the facility must fill out the "Permit to bring animals into the Centre for Cellular Imaging" form and follow the guidelines in that document.
- **Biohazards and chemical hazards**
  - If you work with any hazardous (contagious, carcinogenic, radioactive etc.) material ask the CCI staff if you are allowed to bring these into the facility and, if so, what precautions to take.
  - Any surfaces that have come into contact with specimens that present a biosafety concern must be wiped off with the ethanol solution available in the facility. These surfaces include the instrument itself, the stage, holders, the focus knob, buttons, keyboard as well as the working areas around the instrument.
  - Any waste, which needs to be taken care of as a biohazard (e.g. chemicals, human tissue, infected or genetically modified material), has to be brought back to the user's own lab!
  - Hazardous substances used in the CCI preparation lab have to be handled and disposed of as outlined in the EM lab training, EM lab Safety Handout, MSD sheets and risk assessments.
  - Gloves are prohibited when using any of the microscopes in the CCI facility (**unless temporary additional rules state otherwise!**). All instruments are considered 'clean' and wearing gloves might contaminate them. If gloves are needed for handling specimens they must be removed before operating any instrument. **Handling common equipment with contaminated gloves places all users of the facility at risk of exposure!**
- **Food and drinks are not allowed into any part of the CCI Facility.** If you need a refreshment while working in microscope room or preparation labs, please use one of the common lunch-rooms at the Faculty

**These rules and regulations are subject to change at any time. Suggestions are always welcome. Thank you very much for your cooperation, and do not hesitate to contact the CCI staff in case of doubt!**

Julia Fernandez-Rodriguez  
Head of the Centre for Cellar Imaging